

BE-50 LED

HIGHLIGHTS

- ✓ Professional stereo microscopes for extra long working distance
- ✓ Fixed magnifications from 2.5x up to 20x
- ✓ 0.5x and 1x objectives
- ✓ Overhanging stand
- ✓ LED illumination
- ✓ 2 years warranty



Specifications

Eyepieces

Pair WF 10x/20 mm

Head

Binocular 45° inclined tube

Interpupillary distance adjustable between 55 and 75 mm

Objectives

Standard delivered with 0.5x fixed magnification objective, magnification 5x, field of view 10 mm, working distance 250 mm

Additional 1x fixed magnification objective available, magnification 10x, field of view 20 mm, working distance 230 mm

Other magnifications available with optional eyepieces

All optical parts are anti-fungus treated

Stand

Stable overhanging stand with pillar

Focusing adjustment

Coarse adjustment with tension control

Illumination

Incident 1 W LED illumination built in stand. Intensity is adjustable

Models	Feature				
	Bino	0,5x objective	1,0x objective	Working distance	LED illumination
BE.1802	•	•		250 mm	•
BE.1812	•		•	230 mm	•

Accessories & Spare parts

BE.6005 Pair HWF 5x/22 mm eyepieces
 BE.6010 Pair HWF 10x/20 mm eyepieces
 BE.6015 Pair HWF 15x/12 mm eyepieces
 BE.6020 Pair HWF 20x/10 mm eyepieces

BE.6110 HWF 10x/20 mm eyepiece with micrometer
 BE.6099 Pair eyeshades

BE.7405 0.5x objective for BE-50 LED. Working distance 250 mm
 BE.7410 1x objective for BE-50 LED. Working distance 230 mm

AE.1112 Object micrometer 50 mm divided in 500 parts on glass slide 76 x 26 mm

