

BlueLine
MADE IN HOLLAND

by euromex

ENGLISH



OBSERVE WITH BLUE

 **euromex** *Since 1966!*
microscopes holland

EcoBlue

- Economical microscopes for education
- Monocular models
- Digital and polarisation models available
- LED illumination
- Ergonomical carrying grip
- 5 years warranty



Model EC.1001

The EcoBlue series offer quality biological microscopes at affordable prices. These modern microscopes were specifically developed for education with a special attention for ergonomics and with an outstanding price/quality ratio

FEATURES

Model	Monocular	Digital	S100x obj	S60x obj	LED	Polarisation
EC.1001	•				•	
EC.1101	•		•		•	
EC.1601	•			•	•	
EC.1005	•	•			•	
EC.1105	•	•	•		•	
EC.1605	•	•		•	•	
EC.2001-POL	•					•
EC.2101-POL	•		•			•
EC.2601-POL	•			•		•

Digital



Polarisation



EduBlue

- Economical binocular models
- Dual and Triple magnification objectives
- Available in rack & pinion and pillar stands
- Digital models available
- LED illumination
- Ergonomical carrying grip
- 5 years warranty



Rack & pinion stand
Model ED.1402-S



Pillar stand
Model ED.1402-P

The EduBlue stereo microscopes are specifically designed for educational applications and come with both dual and triple magnifications

FEATURES

Model	Digital head	1x/3x objective	2x/4x objective	1x/2x/3x objective	1x/2x/4x objective	Pillar stand	Rack & Pinion stand
ED.1302-P		•				•	
ED.1302-S		•					•
ED.1402-P			•			•	
ED.1402-S			•				•
ED.1405-S	•		•				•
ED.1502-S				•			•
ED.1802-S					•		•
ED.1805-S	•				•		•

StereoBlue

- Dual 1x/3x or 2x/4x magnification versions
- Zoom Stereo 0.7x to 4.5x versions
- Ergonomically designed stand
- 3 W LED illuminators
- Ergonomical carrying grip
- 5 years warranty



SB.1903



SB.1402

The StereoBlue stereo microscopes are available in fixed dual magnification and zoom 0.7x to 4.5x magnification versions. For education, the amateur and suitable for many routine applications for the industry

FEATURES


Model	Binocular head	Trinocular head	1x/3x objective	2x/4x objective	0.7x-4.5x Zoom
SB.1302	•		•		
SB.1402	•			•	
SB.1902	•				•
SB.1903		•			•

BioBlue

- Microscopes for education
- Monocular and binocular
- Semi-Plan objectives and a reversed quadruple nosepiece
- Integrated mechanical stage
- Digital and polarisation models available
- NeoLED illumination
- Cordless operation
- 10 years warranty



Binocular model BB.4260



Monocular model BB.4200

The beautiful stylish design of the BioBlue series stands out among all other microscopes of this kind. It has been especially designed for biology classes of educational institutes and small laboratories and offer a full range of types

FEATURES

Model	Monocular	Binocular	Digital head	S60x objective	S100x objective	NeoLED	Mechanical stage	Polarisation
BB.4200	•					•		
BB.4220	•					•	•	
BB.4240	•			•		•	•	
BB.4250	•				•	•	•	
BB.4260		•			•	•	•	
BB.4205	•		•			•		
BB.4225	•		•			•	•	
BB.4245	•		•	•		•	•	
BB.4255	•		•		•	•	•	
BB.4267		•	•		•	•	•	
BB.4220-POL (1)	•							•
BB.4240-POL (1)	•			•				•
BB.4260-POL (1)		•			•			•

(1) Polarisation models are standard available as halogen version



by euromex



Euromex Microscopen b.v.

is a leading manufacturer of microscopes and other optical instruments. Founded in 1966, Euromex has become a world-class supplier of biological and stereo microscopes

The corporate office is based in Arnhem, The Netherlands. A facility with a 2.000 m² conditioned logistics warehouse, an opto-mechanical workshop, an R&D department and a high-level quality control department

Around the world, Euromex operates in more than 80 countries through distributors, resellers and agents. A wide variety of customers such as schools and educational institutes, clinical and research laboratories and a broad range of industrial customers are using Euromex microscopes



• ISO 9001

The Euromex Quality System is certified according to ISO 9001:2008 and supports our pursuit of continuous improvement and our on-going commitment to provide our world-wide customers assurance of product quality



Euromex microscopen

info@euromex.com
www.euromex.com

Papenkamp 20
6836 BD Arnhem
Netherlands
Tel: +31 (0) 26 323 22 11

Carretera de Barcelona 88, Entresuelo
Edificio Technomar
08302 Mataró , Spain
Tel: +34 717 168 377